

Analyst®

Software
Version 1.6.3 HotFix 4

Fixed Issue

User-Created Formula columns in the Results Table are not automatically updated when samples are removed from the Results Table

In the Quantitation module of the Analyst® 1.6.3 software and the Analyst® 1.6.3 software with HotFix 3, if one or more samples are removed from the Results Table, then user-created Formula columns in the Results Table are not automatically updated. As a result, the content in the Formula cells becomes out of sync with the sample entry rows that follow the deleted samples.

This incorrect data is presented in the Results Table. The data can be copied using Ctrl + C and then exported to a text or .pdf file, or printed. (AN-1311)

HotFixes and Patches Included in the Analyst® 1.6.3 HotFix 4

- Analyst® 1.6.3 HotFix 1
- Analyst® 1.6.3 HotFix 2
- Analyst® 1.6.3 HotFix 3
- Analyst® 1.6.3 Patch for Negative CE Issue (also included in Analyst® 1.6.3 HotFix 2 and Analyst® 1.6.3 HotFix 3)
- Analyst® 1.6.3 Patch for Reporter Analyte Ratio Plots (also included in Analyst® 1.6.3 HotFix 3)
- Analyst® 1.6.3 Patch for ExionLC (also included in Analyst® 1.6.3 HotFix 3)
- Analyst® 1.6.3 Patch for Audit Trail Issue (also included in Analyst® 1.6.3 HotFix 3)
- Analyst® 1.6.3 Patch for Peak Review (also included in Analyst® 1.6.3 HotFix 3)

Where to Get Help

Other Documentation

- Analyst® 1.6.3 Software Release Notes
- Analyst® 1.6.3 Software Installation Guide
- Analyst® 1.6.3 HotFix 1 *Release Notes*
- Analyst® 1.6.3 HotFix 2 *Release Notes*

- Analyst® 1.6.3 HotFix 3 *Release Notes*

Support

- sciex.com/request-support
- sciex.com/contact-us

Requirements

- Analyst® 1.6.3 software

Note: The Analyst® 1.6.3 HotFix 4 installer silently removes any HotFix and patch listed in [HotFixes and Patches Included in the Analyst® 1.6.3 HotFix 4](#) if it is installed on the computer.

Install the Software

This HotFix is compatible with Windows 7 (32- or 64-bit) operating systems.

1. Log on to the computer as a user with Administrator privileges.
2. Deactivate the hardware profile and then close the Analyst software.
3. Extract the content of **Analyst163HF4.zip** file to your hard drive.
4. Navigate to this location on the hard disk, double-click the **setup.exe** file and then follow the on-screen instructions.

When the HotFix is installed, a new entry is added to the Programs and Features list in Control Panel. Additionally, a shortcut to these release notes is placed under **Start > All Programs > SCIEX > Analyst**.

Updated files

The HotFix adds the *Analyst 1.6.3 HotFix 4 Release Notes.pdf* to the:

- C:\Program Files\Analyst\Help folder (on Windows 7, 32-bit computers).
- C:\Program Files (x86)\Analyst\Help folder (on Windows 7, 64-bit computers).

The HotFix updates the *QuantRT.ocx* file located in the:

- C:\Program Files\Analyst\Bin folder (on Windows 7, 32-bit computers)
- C:\Program Files (x86)\Analyst\Bin folder (on Windows 7, 64-bit computers)

Remove the HotFix

1. Deactivate the hardware profile and then close the Analyst software.
2. Open **Control Panel** and then click **Programs and Features**.
3. Select **Analyst® 1.6.3 HotFix 4** from the list and then click **Uninstall**.

The HotFix is removed from the program list. The *Analyst 1.6.3 HotFix 4 Release Notes.pdf* file and shortcut are also removed. The *QuantRT.ocx* file is reverted to

the original state for the Analyst[®] 1.6.3 software.

Known Issues and Limitations

None

This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate those manufacturers' products as supplied by SCIEX for incorporation into its equipment and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are SCIEX's sole and exclusive representations, warranties, and obligations. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.

For research use only. Not for use in diagnostic procedures.

AB Sciex is doing business as SCIEX.

The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners.

AB SCIEX™ is being used under license.

© 2017 AB Sciex



AB Sciex Pte. Ltd.
Blk 33, #04-06
Marsiling Ind Estate Road 3
Woodlands Central Indus. Estate
SINGAPORE 739256

