



Answers for Science.  
Knowledge for Life.™



### Shane Needham

Ph.D., Co-Founder and CEO,  
Alturas Analytics, Inc.

## Case Study

### Project Goal

To assess the advantages of moving from ligand binding assays (LBA) to microflow LC-MS/MS for large molecule bioanalysis and how the source design affects performance gains.

### The Challenges

- Inadequate sensitivity and selectivity for some bioanalytical assays using LBA or conventional flow LC-MS/MS
- Any new equipment and methodologies must be rugged and reliable, and require minimal optimization
- Any changes to established equipment and methodologies must require minimal new training of lab personnel

### The Solution

- Microflow LC-MS/MS using the Optiflow® Ion Source with our SCIEX triple quadrupole instruments. Column directly connects to source.

### The Outcomes

- Large gains in peak area and peak height with microflow LC-MS/MS and Optiflow Source for large and small molecule bioanalysis, equating to major improvements in S/N (4-20x) and lower LLOQ
- Direct connection port for LC column minimizes dead volume resulting in less band broadening, less tailing, and more narrow peaks
- Source requires minimal to no optimization
- Source requires less than 10 minutes of training
- Have now analyzed many samples to support both clinical and pre-clinical large and small molecule studies by microflow LC-MS/MS using the Optiflow Source including some challenging programs that had selectivity and sensitivity issues that were solved with this workflow

“How a source is designed for microflow is critical for performance. You can’t just use a conventional source at microflow rates and expect to get good results. We love the Optiflow source. You put it on and it works. Optiflow is well-designed, convenient, easy to set up in seconds, and we get better results. No adjustments necessary. That convenience is worth gold.”

### Type of Organization

Contract Research Organization (CRO) specializing in pre-clinical and clinical LC-MS/MS and GC-MS/MS bioanalysis.

### Goals

Development and validation of assays for the determination of therapeutic agents and biomarkers from biological matrices, supporting drug candidates from discovery through phase IV clinical trials.

### Applications

Bioanalysis, method development and validation, and biomarker quantitation.

### SCIEX Products

- Triple Quad Systems including 6500+, 6500, 5500, API 5000, 4000
- Optiflow™ Turbo V Ionization Sources
- Optiflow™ Interfaces

For Research Use Only. Not for use in diagnostic procedures.

AB Sciex is operating as SCIEX. © 2019 AB Sciex. The trademarks mentioned herein are the property of the AB Sciex Pte. Ltd. or their respective owners. AB Sciex™ is being used under license.

RUO-MKT-07-9298-A

Answers for Science. Knowledge for Life.™