

Clinical applications of direct coupling of Bio-compatible SPME devices to MS via Open-Port Probe sampling interface

M. Tascon¹, G. Gomez-Rios¹, C. Liu², A. Roszkowska¹, Nikita Looby¹, N. Reyes-Garcés¹, D. Arnold², T. Covey³, B. Bojko^{1, 4}, J. Pawliszyn¹

¹Department of Chemistry, University of Waterloo, Ontario, Canada, N2L 3G1, ²SCIEX, 71 Four Valley Drive, Concord, Ontario L4K 4V8, Canada, ³SCIEX, 1201 Radio Road, Redwood City, CA 94065, USA, ⁴Department of Pharmacodynamics and Molecular Pharmacology, Faculty of Pharmacy, Nicolaus Copernicus University in Torun, Dr. A. Jurasza 2, 85-089 Bydgoszcz, Poland.


Abstract

We are presenting the advances achieved up to date employing a microfluidic open port that allows the coupling of biocompatible solid-phase microextraction (Bio-SPME) devices through the generation of a stagnant volume followed by a three second switch flow. In this way, an efficient desorption and very sharp peaks (2-3 seconds FWHM) are reached.


Towards the monitoring of drugs and biomarkers in the operation room, some applications were developed using this interface such as the monitoring of doxorubicin from perfusate used during in-vivo (IVLP) and ex-vivo (EVLP) lung perfusion to treat lung cancer.

Regarding the high-throughput therapeutic drug monitoring, the aim of this work is to push down the limits of detection trying to improve dramatically the total analysis time in high-throughput and the individual turn-around times, as well. Therefore, advances in the development of the quantitation method of testosterone from human plasma is herein presented. The concept of rapid on-fiber derivatization is gathered by using the Amplifex™ keto reagent. Finally, quantitation of immunosuppressive drugs from 100 µL of whole blood is presented. LOQs of sub-ppb levels were achieved for all the compounds with turn-around times of less than 90 minutes.


Instrumental




Set-up and Fundamentals



Analytical workflow




Instrumental cycle




Results


Monitoring of doxorubicin from ex-vivo perfusate from EVLP and IVLP




Quantitation of Testosterone from human plasma



On-fiber derivatization



Quantitation of immunosuppressive drugs from whole blood



Conclusions

- The new microfluidic open port herein presented showed to be a highly sensitive interface for direct coupling of SPME to MS via ESI ionization.
- Doxorubicin was successfully quantitated from real perfusate samples from different surgical procedures obtaining a LOQ of 1 ppm.
- A fast and sensitive method for the determination of immunosuppressive drugs from 100 µL of whole blood was developed. Turn-around times of individual samples is less than 90 min.
- Testosterone was directly quantitated from human plasma obtaining 10 ppb as LOQ. In order to improve this performance, an on-fiber derivatization procedure is presented showing off-line preliminary results with an LOQ close to 1 ppb.

References

- Gomez-rios, G. A.; Liu, C.; Tascon, M.; Reyes-Garcés, N.; Arnold, D. W.; Covey, T. R.; Pawliszyn, J. Anal. Chem. 2017, 89 (7), 3805–3809.
- Star-Weinstock, M.; Williamson, B. L.; Dey, S.; Pillai, S.; Purkayastha, S. Anal. Chem. 2012, 84 (21), 9310–9317.

Acknowledgments