



# Analysis of Nitrofurans Metabolites in Honey Using the SCIEX Triple Quad™ 3500 System

Chandrasekar<sup>1</sup>. M; Santosh Kapil.G<sup>1</sup>; Aman Sharma<sup>1</sup>; Anoop Kumar<sup>1</sup>; Manoj G Pillai<sup>1</sup> & Jianru Stahl-Zeng<sup>2</sup>  
<sup>1</sup>SCIEX, Udyog Vihar, Phase IV, Gurgaon - 122015, Haryana, India; <sup>2</sup>SCIEX, Darmstadt, Germany

## Overview

A LC-MS/MS method for the simultaneous quantification of four Nitrofurans Metabolites (3-amino-2-oxazolidinone (AOZ) for furazolidone, 3-amino-5-methylmorpholino-2-oxazolidinone (AMOZ) for furaltadone, 1-aminohydantoin (AHD) for nitrofurantoin and semicarbazide (SEM) for nitrofurazone) on SCIEX Triple Quad™ 3500 was developed to detect Nitrofurans residues in honey samples. The method showed adequate linearity with correlation coefficients above  $r \geq 0.99$  for all four analytes. The Minimum Required Performance Limit (MRPL) for Nitrofurans Metabolites in Honey was  $1 \mu\text{g}/\text{kg}$ .

## Introduction:

Nitrofurans are broad spectrum antibacterial agents which were used in the treatment of bacterial infections in bee colony health. Nitrofurans have been prohibited in food-producing animals in the European Union and most other Countries for public health and safety concerns. The nitrofurans are unstable and easily metabolized within a few hours but Nitrofurans metabolites are highly stable in nature. Several methods have been described in the analysis of nitrofurans metabolite in honey samples by incubation period for derivatization with nitrobenzaldehyde in overnight or 16 hours at  $37^\circ\text{C}$ .

The LC-MS/MS method developed on SCIEX Triple Quad™ 3500 described here for the quantitation of nitrofurans metabolites in honey was found to meet the regulatory requirements of  $1 \mu\text{g}/\text{kg}$ .



Figure 1: SCIEX Triple Quad™ 3500



## LC Conditions

LC separation was achieved using the Shimadzu prominence system with an Eclipse plus C18 (4.6x150 mm) 5 µm column with a gradient of 1mM ammonium acetate as mobile phase A and Methanol as mobile phase B at flow rate of 0.4 mL/min. The injection volume was set to 10 µL. Gradient profile is given Table1.

Time (min)	Mobile phase A%	Mobile phase B%
0.01	95	5
0.50	45	55
3.50	45	55
4.00	95	5
12.00	Controller	Stop

Table 1: Mobile Phase Gradient

## MS/MS Conditions

The SCIEX Triple Quad™ 3500 was operated in Multiple Reaction Monitoring (MRM) mode. The Turbo V™ source was used with an Electrospray Ionization (ESI) probe in positive polarity. Two selective MRM transitions were monitored for all nitrofurantoin metabolites using the Analyst® 1.6.2 Software and MultiQuant™ Software version 3.0.2. MRM transition is given in Table 2.

Compound	Precursor ion	Product ion Quantifier	Product ion Qualifier
AOZ	236.0	104.0	78.0
AMZO	335.0	291.1	128.2
SEM	209.0	166.0	192.0
AHD	249.1	134.0	104.0

Table 2: MRM transitions

## Results and Discussions

The results of repeatability data obtained for nitrofurantoin metabolites in the honey matrix is given in table 3 at different levels.




Figure 2: Structures of Nitrofurantoin, Nitrofurantoin metabolites and Nitrophenyl derivatives.

## Materials and Methods

### Chemicals

Nitrofurantoin metabolite standards were purchased from clearsynth and 2-Nitrobenzaldehyde was purchased from Sigma Aldrich ≥99% Purity. All other chemicals used were of LC-MS grade, commercially available.

### Honey samples

Honey samples were procured from local market of Delhi and Gurgaon, India and were stored at room temperature until end of analysis.

### Sample Preparation

Honey sample (1gm) was mixed with 3ml of HCl (0.1M) and 50mM of 2-Nitrobenzaldehyde (0.3ml), vortexed and incubated on ultrasonic bath for 16hr added 0.6ml of 1M K2HPO4 solution and added 10 ml of ethyl acetate, vortexed it, followed by centrifugation at 4000 rpm. The supernatant was evaporated to dryness, reconstituted with 1ml of Methanol: water (5:95) and 10µl is used for LC-MS/MS analysis.




Figure 3: Signal to noise of AMZO, AHD, SEM and AOZ at MRPL level (1.0 ng/ml) in honey matrix sample

For all four Nitrofurantol metabolites in honey the matrix based calibration curve shows excellent linearity (0.50 to 20.0 ppb), with a correlation coefficient  $r \geq 0.99$  using linear regression and weighing factor  $1/x$ . The SCIEX Triple Quad™ 3500 was found to be capable of analyzing concentrations well below the MRPL required by EU. The signal to noise ratio for all four nitrofurantol metabolites at 1.0 ppb is  $\geq 30$ . The signal to noise ratios and calibration curves are shown in Figure 3 and Figure 4.

Analyte	Repeatability			Recovery (n=6)		
	½ MRPL (0.5ppb)	MRPL (1.0ppb)	1.5MRPL (1.5ppb)	½ MRPL (0.5ppb)	MRPL (1.0ppb)	1.5MRPL (1.5ppb)
AOZ	6.01	7.00	4.28	113.47	95.05	89.89
AMZO	12.16	4.46	4.33	83.80	103.88	96.11
SEM	4.49	7.31	9.44	109.67	98.13	91.33
AHD	4.96	7.58	8.22	114.90	105.40	105.00

Table 3: Repeatability (%CV) and recovery statistics in honey sample

Analyte	Calibration Range (ppb)	Linearity (r)	CC $\alpha$	CC $\beta$
AOZ	0.5 -20	0.9981	0.58	0.63
AMZO	0.5 -20	0.9963	0.62	0.70
SEM	0.5 -20	0.9974	0.56	0.60
AHD	0.5 -20	0.9987	0.57	0.16

Table 4: Summary of CC $\alpha$ , CC $\beta$  and linearity in honey sample

Recovery experiments were performed in honey samples at ½ MRPL, MRPL and 1.5 MRPL level (n=6). The recovery of all nitrofurantol metabolites was  $\geq 80\%$ . The recovery data for nitrofurantol metabolites are shown in Table 3. The retention time (RT) of the AHD, AOZ, SEM, AMZO, were 5.94, 6.57, 6.57 and 7.50 min, respectively.




Figure 4: Matrix based calibration curve AOZ, AMZO, SEM and AHD in Honey sample showing  $r = >0.99$

Decision limit ( $CC\alpha$ ) and detection capability ( $CC\beta$ ) were calculated for all the four derivatives of Nitrofurans in Honey samples. The calculation was based on using linear regression model analyzing spiked honey samples at below MRPL level (Van Loco et al, 2007).

The calculated value of  $CC\alpha$  and  $CC\beta$  are given in table 4. The decision limit ( $CC\alpha$ ) and detection capability ( $CC\beta$ ) of all the metabolites were well below the MRPL.

## Conclusions

- The developed quantitative method of Nitrofurans in honey on SCIEX Triple Quad™ 3500 was sensitive, linear, and reproducible.
- Trueness (Average recovery %) for this method found to be  $\geq 80\%$  at various MRPL levels.
- The method and data presented in the application note showcase the fast and accurate solution for the quantitation and identification of nitrofurans metabolites in honey samples for quality control check

## References

1. U.S. Food and Drug Administration Center for Food Safety Applied Nutrition Food Compliance Program Chapter 03 – Foodborne Biological Hazards (10-01-97)  
<http://www.cfsan.fda.gov/~comm/cp03039.html> accessed 2/13/09
2. Fatih Alkan, Arzu Kotan, Nurullah Ozdemir; Development and validation of confirmatory method for analysis of nitrofurans metabolites in milk, honey, poultry meat and fish by liquid chromatography-mass spectrometry. Macedonian Veterinary Review 39 (1): 15-22 (2016)
3. Nurullah Ozdemir, Fatih Alkan, Arzu Kotan; Rapid confirmatory method for analysis of nitrofurans metabolites in egg by liquid chromatography-mass spectrometry. International Journal of Technical Research and Applications 4(2): 31-37 (2016)
4. <http://www.eicindia.gov.in/Services/Pre-Compliance/Residue-Monitoring-Plans.aspx>
5. EC (2010): Council Regulation 37/2010/EU of 22 December 2009 on pharmacologically active substances and their classification regarding maximum residue limits in foodstuffs of animal origin. Official Journal of European Union, L15, 1-72.
6. Tomasz Śniegocki, Andrzej Posyniak, And Jan Żmudzki, Determination of Nitrofurans Metabolite Residues In Eggs By Liquid Chromatography-Mass Spectrometry. Bull Vet Inst Pulawy 52, 421-425 (2008)
7. J. Van Loco, A. Janosi, S. Impens, S. Fraselle, V. Cornet, J.M. Degroot, Calculation of the decision limit ( $CC\alpha$ ) and the detection capability ( $CC\beta$ ) for banned substances: The imperfect marriage between the quantitative and the qualitative criteria. Analytica Chimica Acta Volume 586, Issues 1–2, 14, Pages 8–12 (2007)

AB Sciex is doing business as SCIEX.

© 2017 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT-02-5033-A