

LC Solutions Tailored to Meet Your Needs

Complete. Reliable. Expandable.



Answers for Science.
Knowledge for Life.™

Your LC-MS laboratory has unique analytical requirements. SCIEX provides LC-MS solutions tailored to meet your needs. Introducing the ExionLC™ Series of robust and reliable HPLC/UHPLC systems offering a wide range of performance and flexibility.



Expandability, Robustness
and Ultra-low Carryover



Reliable and Expandable
Platform with High Uptime



ExionLC™ AD

The ExionLC AD UHPLC offers expandability, robustness and ultra-low carryover performance for the most challenging method development and sample analysis applications – covered by 40 years of trusted LC-MS innovation and SCIEX service.

- Near zero carryover with 3 different rinse solutions plus separate rinse port solution
- Injection cycle time as low as 14 seconds per sample
- Sample cooling, 2 microtiter plates or 105 vials
- Pressures to 18,850 PSI (1,300 bar)
- Flow rates from 0.0001 to 10.0 mL/min
- Heated block oven to 150°C
- Optional accessories include 6 well plate autosampler, rack changer, and optical detectors


ExionLC™ AC

A reliable and expandable HPLC/UHPLC system with the power of industry-proven SCIEX OS or Analyst® Software, the ExionLC™ AC ensures very short turn-around time to results with ultra-low downtime across thousands of injections. By combining streamlined quality methods, industry-proven software, and leading edge LC-MS/MS technology, you can now quickly take advantage of ultimate sample throughput, accuracy, and reliability.

- Modular pumping system
Up to 3 mL/min. at 9,500 PSI (650 bar)
Up to 5 mL/min. at 6,000 PSI (400 bar)
- Sample cooling, 2 microtiter plates or 105 vials
- Large capacity solvent tray easily holds up to 2 liter bottles for uninterrupted analyses
- Solvent degasser with 5-channels
- Forced air column oven with temperature control from 10°C below ambient to +85°C
- Options – rack changer, optical detectors

Network Your Way to Success

Discover How the SCIE X Success Network Speeds and Simplifies Your Path to Answers



Start Your Path to Success Now: [sciex.com/support](https://www.sciex.com/support)

Answers for Science.
Knowledge for Life.™

AB Sciex is doing business as SCIE X.

© 2019 AB Sciex. For Research Use Only. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIE X™ is being used under license.

RUO-MKT-03-2524-B 01/2019

Headquarters

500 Old Connecticut Path
Framingham, MA 01701 USA
Phone 508-383-7700
[sciex.com](https://www.sciex.com)

International Sales

For our office locations please call the division headquarters or refer to our website at [sciex.com/offices](https://www.sciex.com/offices)

