

## 2 天用户现场高级 LC-MS 定量与故障排查培训所需耗材

SCIEX 培训课程遵循经验证的间隔学习法，以最大限度地提高学习效果和记忆。该培训过程在用户现场进行，是讲师指导、实验室上机实操和自主在线学习的独特结合。本培训提供一系列故障排查练习、操作规范，最后以一个对学习者的问题解决测验作为结尾。

为了提升培训效果，需要在培训前安装和配置好 LC-MS 系统。表 1 列出了培训中所需用到的耗材。

表 1: 培训所需耗材

描述	产品编码 (Part number)	规格
三嗪标准溶液 (Triazine standard solution)	4376887	N/A
飞诺美 (Phenomenex) Synergi 4 $\mu$ m Fusion-RP 80Å HPLC 色谱柱	4376878	50 mm $\times$ 2.0 mm
PEEK 管	4425163	0.13 $\times$ 3000 mm
PEEK 管切刀	011281	N/A
注射器	WC010615	1 mL
注射器针头	1005819	1 mL
注射器转接头	1008236	N/A
LC/MS 标准品套装 (含 High/Low PPG 溶液)	4406127	N/A
PEEK 三通接头 (Fitting PEEK Tee 0.020 in BORE)	1006550	N/A
LC-MS 级水	N/A	N/A
MS 级甲醇 (保存于玻璃瓶)	N/A	N/A
MS 级乙腈 (保存于玻璃瓶)	N/A	N/A
LC-MS 级异丙醇	N/A	N/A
甲酸	N/A	N/A
无尘纸	018027	11 x 21 cm, 4.3 x 8.3 英寸
小聚酯棉签	1017396	N/A
Alconox 清洁包 (Alconox packets)	400667DHR	N/A

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to [www.sciex.com/diagnostics](http://www.sciex.com/diagnostics). All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries (see [www.sciex.com/trademarks](http://www.sciex.com/trademarks)).

© 2021 DH Tech. Dev. Pte. Ltd. GEN-CST-02-7730-ZH-C



**Headquarters**  
500 Old Connecticut Path | Framingham, MA 01701 USA  
Phone 508-383-7700  
[sciex.com](http://sciex.com)

**International Sales**  
For our office locations please call the division  
headquarters or refer to our website at  
[sciex.com/offices](http://sciex.com/offices)